

M3-D Agit. ARGOLab

magn.heater square ceramic

Code: 22003013



TECHNICAL DATA SHEET

Performance

Heating power	1000
Over Temperature Safety Device	580
Working temperature	Ta ... 550
Heating Temperature Range	Ta ... 550 °C
Working humidity	80 % RH
Stirring speed	0 ... 1500
Max stirring quantity	10 L

Features

IP rating	21
Temperature indication	LED Display
Speed indication	Digital
Maximum length of magnetic anchor	80
Working plate - material	Glass ceramic

Feeds

Power supply	220
Frequency	50 / 60

Dimensions and weight

Dimensions (WxDxH)	215 x 360 x 112
Weight	4.5
Working plate - dimensions	Ø135 mm