

ArgoLab combined fridge-freezer FRZ 300G - glass door

Code: 41500112



TECHNICAL DATA SHEET

Performance

temperature measurement range	2 ... 8 / -25 ... -10 °C
Recovery time 5°C	19 m
Temperature variation at 5 °C	1,2 °C
Gradient at 5 °C	1.3 °C
Recovery Time at -20°C	26 m
Temperature variation at -20 °C	1 °C
Gradient at -20°C	1.8 °C

Features

Door type	Blind + Glass
Maximum shelf load	15 kg
Through-hole diameter	10 mm
Display	microprocessor/display
Defrosting function	Automatic (FR) and manual (FZ)
Internal light	Yes (FR)
Number of shelves / baskets (only FZ 300-P)	5
Refrigeration	Dynamics (FR) and Static (FZ)
Useful volume	300 L

Feeds

Power supply	220 - 240 V
Power	340 W
Energy consumption per year	548 KW/year
Energy consumption measurement method	Measurement at empty load +5°C / Ambient temperature +25 °C

Dimensions and weight

External dimensions (WxDxH)	600 x 610 x 1770 mm
Internal dimensions (WxDxH)	FR: 510 x 490 x 730 / FZ: 450 x 460 x 440 mm

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

Packaging dimensions (WxDxH)	850 x 790 x 1970 mm
Noisiness	55 dB
Minimum distance between shelves	13 mm
Weight	75 kg
