

ArgoLab Refrigerator FR 110 - solid door

Code: 41500012



TECHNICAL DATA SHEET

Performance

temperature measurement range	2 ... 8 °C
Recovery time 5°C	12 m
Temperature variation at 5 °C	1,3 °C
Gradient at 5 °C	1.3 °C

Features

Door type	Blind
Maximum shelf load	15 kg
Through-hole diameter	10 mm
Display	microprocessor/display
Defrosting function	Automatic and manual
Internal light	Yes
Number of shelves / baskets (only FZ 300-P)	3
Refrigeration	Dynamics
Useful volume	110 L

Feeds

Power supply	220 - 240 V
Power	120 W
Energy Consumption per 24 hours	1.5 kWh/24h
Energy consumption per year	584 KW/year
Energy consumption measurement method	Measurement at empty load +5°C / Ambient temperature +25 °C

Dimensions and weight

External dimensions (WxDxH)	600 x 560 x 805 mm
Internal dimensions (WxDxH)	500 x 420 x 570 mm
Noisiness	55 dB
Minimum distance between shelves	13 mm

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

Weight 50 kg
