

ArgoLab Refrigerator FR 260 G

- glass door

Code: 41500042



TECHNICAL DATA SHEET

Performance

temperature measurement range	2 ... 8 °C
Recovery time 5°C	18 m
Temperature variation at 5 °C	2,2 °C
Gradient at 5 °C	2.6 °C

Features

Door type	Glass
Maximum shelf load	15 kg
Through-hole diameter	10 mm
Display	microprocessor/display
Defrosting function	Automatic and manual
Internal light	Yes
Number of shelves / baskets (only FZ 300-P)	6
Refrigeration	Dynamics
Useful volume	260 L

Feeds

Power supply	220 - 240 V
Power	160 W
Energy Consumption per 24 hours	2.1 kWh/24h
Energy consumption per year	548 KW/year
Energy consumption measurement method	Measurement at empty load +5°C / Ambient temperature +25 °C

Dimensions and weight

External dimensions (WxDxH)	600 x 560 x 1525 mm
Internal dimensions (WxDxH)	500 x 450 x 1290 mm
Noisiness	55 dB
Minimum distance between shelves	13 mm

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

Weight 75 kg
