

ArgoLab Freezer FZ 260

Code: 41500232



TECHNICAL DATA SHEET

Performance

temperature measurement range	-25 ... -10 °C
Recovery Time at -20°C	30 m
Temperature variation at -20 °C	2.1 °C
Gradient at -20°C	2.1 °C

Features

Door type	Blind
Maximum shelf load	15 kg
Through-hole diameter	10 mm
Display	microprocessor/display
Defrosting function	Manual with dedicated key
Internal light	Yes
Number of shelves / baskets (only FZ 300-P)	6
Refrigeration	Statics
Useful volume	260 L

Feeds

Power supply	220 - 240 V
Power	190 W
Energy Consumption per 24 hours	2.3 kWh/24h
Energy consumption per year	949 KW/year
Energy consumption measurement method	Measurement at empty load -20°C / Ambient temperature +25 °C

Dimensions and weight

External dimensions (WxDxH)	600 x 560 x 1525 mm
Internal dimensions (WxDxH)	500 x 450 x 1290 mm
Noisiness	55 dB
Minimum distance between shelves	13 mm

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

Weight 75 kg
