

ArgoLab Chest Freezer FZ 300-P

Code: 41500302



TECHNICAL DATA SHEET

Performance

temperature measurement range	-25 ... -10 °C
Recovery Time at -20°C	30 m
Temperature variation at -20 °C	1 °C
Gradient at -20°C	2 °C

Features

Door type	Blind
Maximum shelf load	15 kg
Through-hole diameter	10 mm
Display	microprocessor/display
Defrosting function	Manual with dedicated key
Internal light	No
Number of shelves / baskets (only FZ 300-P)	2
Max number of shelves / baskets (only FZ 300-P)	4
Refrigeration	Statics
Useful volume	285 L

Feeds

Power supply	220 - 240 V
Power	300 W
Energy consumption per year	949 KW/year
Energy consumption measurement method	Measurement at empty load -20°C / Ambient temperature +25 °C

Dimensions and weight

External dimensions (WxDxH)	1280 x 770 x 890 mm
Internal dimensions (WxDxH)	1080 x 480 x 590 mm
Packaging dimensions (WxDxH)	1350 x 810 x 1070 mm
Noisiness	55 dB

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

Weight 69 kg
