

NEYA 10R PROFESSIONAL REFRIGERATED Centrifuge, 4x175 ml, 6000 rpm, 10 Programmes

Code: 40100332



TECHNICAL DATA SHEET

Standards applied

Warranty	2 years
----------	---------

Performance

Capacity max. (fixed angle rotor)	6 x 100 ml
Capacity max. (swing rotor)	4 x 175 ml
Acceleration ramps	0 ÷ 9 (0 = min. - 9 = max.)
Deceleration ramps	0 ÷ 9 (0 = no brake, 1 = min. - 9 = max.)
RCF max (fixed angle rotor)	4800 xg
RCF max (swing rotor)	3600 xg
Noise	55 dB
Maximum speed with fixed angle rotor	6000 rpm
Maximum speed with swing rotor	4500 rpm

Features

Imbalance control	Yes
Time and date	Yes
RCF display	Yes
Spin function	Yes
RCF setting	Yes, from 50xg to 4800xg with steps of 50xg
Rotor indication	Yes
Speed limiter	Yes
Programs	Yes
Control system	Microprocessor and backlit display
Timer	00:30 ÷ 99:50 (mm:ss) e in continuo hh:mm
Type of motor	Brushless induction

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

Max speed check	Yes
-----------------	-----

Feeds

Power supply	220V± 10% - 50 Hz
Consumption	1000 W
Max. amperes	4 A

Dimensions and weight

Entrance height	275 mm
Dimensions (WxDxH)	730 x 640 x 330 mm
Dimensions height when open	670 mm
Weight	70 kg
Weight (with packaging, no rotor)	92 kg