

KERN DLB-A01N.

KERN

Temperature calibration set consists of measuring sensor and display device



Pictograms

STANDARD



Category

Brand	KERN
Product category	Balance
Product group	Temperature calibration set

Measuring System

Units	°C
Temperature [Min] (°C)	40 °C
Temperature [Max] (°C)	200 °C
Readability temperature [d] (°C)	1 °C

Construction

Dimension (W×D×H)	100×0×30 mm
Material housing	metal, lacquered
Cable length	0 m

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	215×200×60 mm
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,30 kg
Shipping weight	0,5 kg

Product information

GTIN/EAN number	4045761188020
Trade number	90160090
REACH classification	no prohibited substances