

KERN DIS-A01

KERN

Temperature calibration set consists of measuring sensor and display device



Environmental conditions

Storage temperature [Min] -10 °C

Storage temperature [Max] 50 °C

Approval

CE mark ✓

Packing & Shipping

Delivery time 1 d

Dimensions packaging (W×D×H) 230×180×45 mm

Shipping method Parcel service

Net weight approx. 0,15 kg

Gross weight approx. 0,50 kg

Shipping weight 0,47 kg

Product information

GTIN/EAN number 404576 1493001

Trade number 90160090

Pictograms

STANDARD



Category

Brand	KERN
Product category	Balance
Product group	Temperature calibration set
Product family	DIS

Measuring System

Units	°C
Temperature [Min] (°C)	-70 °C
Temperature [Max] (°C)	400 °C
Default unit	°C
Readability temperature [d] (°C)	0,1 °C

Construction

Material	plastic + metal
Material housing	ABS plastic
Cable length	1 m

Functions

Auto-off interval(s) in battery mode/rechargeable battery mode	10 min
--	--------

Power Supply

Supplied power supply	Battery
Battery / accumulator type	Alkali(-Manganese)
Battery	3× LR44
Battery connection	Pad + Spring
Battery operating time	80 h
Battery capacity	1.200 mAh
Battery voltage	4,5 V
Battery - Construction	Individual cells
Rechargeable battery operating time - backlight off	80 h