



FRONTIER™ 5515&5515R

High-Speed Micro Centrifuges



Smooth and Dependable High-Speed Micro Centrifuges for Liquid Separation Applications

Frontier 5515 and 5515R Benchtop Micro Centrifuges provide high-speed operation for a wide variety of tube volumes and capacities. With refrigeration capabilities (FC5515R), 10 acceleration/deceleration settings and multiple safety features, the FC5515 and FC5515R ensure stable liquid separation for applications requiring high-speed operation in molecular and cell biology, clinical, research and industrial applications.

Features:

- Excellent functionality when used with rotors compatible with high throughput and BIOSEALS options for 1,5/2,0 ml or 5 ml microtubes, spin columns, and cryo tubes. Additionally, when equipped with the appropriate rotor options, these centrifuges can also be used in diagnostic settings for hematocrit capillaries and PCR tubes/strips.
- Designed for ease of operation, Frontier Micro centrifuges feature convenient control on an intuitive interface, which allows for responsive operation – even when wearing gloves.
- Enclosing all their functionalities within a small footprint, Frontier Micro centrifuges save valuable space on the lab bench.

FRONTIER™ 5515&5515R High-Speed Micro Centrifuges

Easy-to-Use Control Panel

- Easy-to-turn rotary knob for parameter adjustments
- Waterproof keypad provides protection that helps significantly increase the lifespan of the product
- One-touch toggling between speed, G-force, and temperature controls (R model only)
- LCD display with light blue backlight which projects information in real-time
- View and adjust user parameters from the main screen; no need to navigate through sub-screens



High-Performance Driven by Sophisticated Technology

- Powerful induction motor reaches speeds up to 15 200 RPM in less than 16 seconds
- Maximum G-force of 21 953 RCF
- 10 acceleration/deceleration settings
- 99 memory sets to quickly store and recall commonly used parameters
- FC5515R model has refrigeration capabilities from 40 °C to -20 °C
 - Cools down from 23 °C to 4 °C in less than 10 minutes
 - R290, an environmentally friendly coolant, contributes to the eco-friendliness of refrigerated centrifuges by minimizing their environmental impact and supporting efforts to combat global warming
 - Pre-cool feature quickly lowers the temperature of rotor and chamber before a run



Safety Features that Prevent Injury

- Motorized lid-locking system prevents centrifuge from opening during operation
- Rotor imbalance sensing systems stops the unit when it detects an imbalanced load
- Sealed and TÜV certified rotors to safeguard operators from harmful aerosols



Flexibility for Varied Centrifugation Needs

- Works with microtube rotors with different rotor capacity – 24 × 1,5/2,0 ml, 30 × 1,5/2,0 ml or 44 × 1,5/2,0 ml
- Compatible with 6 rotors (7 rotors on FC5515 model), including spin column rotor, hematocrit rotor, PCR Strips/tubes rotor and angle rotor for 5ml tubes
- An ergonomically designed rotors for microtubes, spin columns and PCR Stripes/tubes quickly and easily closes with 1/4 of a turn




FRONTIER™ 5515&5515R High-Speed Micro Centrifuges


Specifications


Model	FC5515		FC5515R	
Order No.	230V: 30130866	120V: 30130867	230V: 83042057	120V: 83042058
Speed Range	200 rpm – 15 200 rpm; 10 rpm / increments			
Maximum RCF	21 953×g; 10×g / increments			
Maximum Capacity (Rotor)	Angle Rotor: 44×1,5/2,0 ml			
Temperature Range	Air cooled		-20 °C to 40 °C, 1 °C / increments	
Refrigerant	N/A		R290 (Propane)	
Running Time	Quick Spin, 10 sec to 99 hr 59 min, or continuous			
Noise Level (depending on the rotor)	≤ 60 + 2 dB(A)		≤ 57 + 2 dB(A)	
Acc/Dec Levels	10/10			
Allowable Density at Maximum Speed	1,2 g/ml			
Allowable Kinetic Energy	7 204 Nm			
Rated voltage	230 V/50-60 Hz	120 V/50-60 Hz	230 V/50-60 Hz	120 V/50-60 Hz
Voltage Fluctuation	± 10%			
Rated current	1,3 A	2,3 A	2,4 A	5,1 A
Power continues	270 W	280 W	500 W	560 W
Dimensions (W×D×H)	280×393×287 mm		280×555×287 mm	
Net Weight (without rotor)	16 kg		35 kg	
Shipping Dimensions (W×D×H)	410×540×400 mm		390×690×390 mm	
Shipping Weight (without rotor)	20 kg		44 kg	
Ambient conditions (EN/IEC 61010-1)	Environment	For indoor use only		
	Altitude	Use up to an altitude of 2000 m		
	Ambient temperature	5 °C up to 35 °C		
	Max. relative humidity	Max. relative humidity 80 % for temperatures up to 31 °C, Decreasing linearly to 50 % relative humidity up to 35 °C.		
EMC	EN/IEC 61326-1 Class B emissions, Basic immunity FCC Class B emissions			
Not suitable for use in hazardous environments.				


Accessories

	Order No.	30372717			
	Description	Rotor Angle 44×1,5/2,0 ml V1			
	Compatible with	FC5515		FC5515R	
	Max. Speed	15 000 rpm		15 000 rpm	
	Max. RCF	21 127×g / 21 379×g		21 127×g / 21 379×g	
Accessories					
Order No.	Description	Vol./Tube	Diameter	Tube type	No. of Tubes/Rotor
Without Adapter	Rotor Angle 44×1,5/2,0 ml V1	1,5/2,0 ml	11 mm	Microtube	44
30130885	Adapter 1×0,5 ml D8 mm, 6 pk	0,5 ml	8 mm	Microtube	44
30130884	Adapter 1×0,2/0,4 ml D6 mm, 6 pk	0,2/0,4 ml	6 mm	Microtube	44


FRONTIER™ 5515&5515R High-Speed Micro Centrifuges


	Order No.	30130872			
	Description	Rotor Angle 30 × 1,5/2,0 ml			
	Compatible with	FC5515		FC5515R	
	Max. Speed	13 500 rpm		13 500 rpm	
	Max. RCF	19 151 × g		19 151 × g	
Accessories					
Order No.	Description	Vol./Tube	Diameter	Tube type	No.of Tubes/Rotor
Without Adapter	Rotor Angle 30 × 1,5/2,0 ml	1,5/2,0 ml	11 mm	Microtube	30
30130885	Adapter 1 × 0,5 ml D8 mm, 6 pk	0,5 ml	8 mm	Microtube	30
30130884	Adapter 1 × 0,2/0,4 ml D6 mm, 6 pk	0,2/0,4 ml	6 mm	Microtube	30


	Order No.	30642371			
	Description	Rotor Angle 24 × 1,5/2,0 ml BIOSEALS V2			
	Compatible with	FC5515		FC5515R	
	Max. Speed	15 200 rpm		15 200 rpm	
	Max. RCF	21 953 × g		21 953 × g	
Accessories					
Order No.	Description	Vol./Tube	Diameter	Tube type	No.of Tubes/Rotor
Without Adapter	Rotor Angle 24 × 1,5/2,0 ml BIOSEALS V2	1,5/2,0 ml	11 mm	Microtube	24
30130885	Adapter 1 × 0,5 ml D8 mm, 6 pk	0,5 ml	8 mm	Microtube	24
30130884	Adapter 1 × 0,2/0,4 ml D6 mm, 6 pk	0,2/0,4 ml	6 mm	Microtube	24

	Order No.	30130873			
	Description	Rotor Angle 12 × 5 ml FA			
	Compatible with	FC5515		FC5515R	
	Max. Speed	14 500 rpm		14 500 rpm	
	Max. RCF	19 978 × g		19 978 × g	
Accessories					
Order No.	Description	Vol./Tube	Diameter	Tube type	No.of Tubes/Rotor
Without Adapter	Rotor Angle 12 × 5 ml FA	5 ml	17 mm	Microtube	12
30130886	Adapter 1 × 1,5/2,0 ml D11 mm, 6 pk	1,5/2,0 ml	11 mm	Microtube	12
30130887	Adapter 1 × 1 ml D12,6 mm Cryo, 6 pk	1 ml	12,6 mm	Cryo Tube - Flat Bottom	12
30130888	Adapter 1 × 1,8 ml D12,6 mm Cryo, 6 pk	1,8 ml	12,6 mm	Cryo Tube - Round Bottom	12

FRONTIER™ 5515&5515R High-Speed Micro Centrifuges

	Order No.	30130874			
	Description	Rotor Angle 4×8-Place PCR Strip			
	Compatible with	FC5515	FC5515R		
	Max. Speed	15 000 rpm	15 000 rpm		
	Max. RCF	15 343×g	15 343×g		
Accessories					
Order No.	Description	Vol./Tube	Diameter	Tube type	No.of Tubes/Rotor
Without Adapter	Rotor Angle 4×8-Place PCR Strip	0,2/0,4 ml	-	8-w PCR Strip	32

	Order No.	83041007			
	Description	Rotor Hematocrit 24×Capillaries V2			
	Compatible with	FC5515	FC5515R		
	Max. Speed	13 000 rpm	-		
	Max. RCF	16 058×g	-		
Accessories					
Order No.	Description	Vol./Tube	Diameter	Tube type	No.of Tubes/Rotor
Without Adapter	Rotor Hematocrit 24×Capillaries V2	-	1,9 mm	Capillaries	24

	Order No.	83041050			
	Description	Rotor Angle 24×2,0 ml Spin Column V1			
	Compatible with	FC5515	FC5515R		
	Max. Speed	13 500 rpm	13 500 rpm		
	Max. RCF	17 317×g	17 317×g		
Accessories					
Order No.	Description	Vol./Tube	Diameter	Tube type	No.of tubes/rotor
Without Adapter	Rotor Angle 24×2,0 ml Spin Column V1	2 ml	11	Spin Column	24

Shortcut list:

ID – auto rotor identification
FA – Conical bottom / Falcon®

Hi – version with higher speed
FB – Flat bottom

Vac – Vacutainer®
RB – Round bottom

Sar – Sarstedt
Rim – Skirted conical bottom

FRONTIER™ 5515&5515R High-Speed Micro Centrifuges



Other Standard Features and Equipment

Microprocessor control, Imbalance switch-off, Audible signal at the end of each run with melody options, Emergency release, 99 memory sets, Ten acceleration and ten deceleration rates, Pre-selection of running time from 10 sec to 99 hr 59 min or continuous, Quick spin for short runs.

Compliance

See Declaration of Conformity available CE Declarations | OHAUS

Outline Dimensions



OHAUS Europe GmbH

e-mail: ssc@ohaus.com

Tel: 0041 22 567 53 19

e-mail: tsc@ohaus.com

Tel: 0041 22 567 53 20

www.ohaus.com

The management system governing the manufacture of this product is ISO 9001:2015 certified.

80774801_E 20260226 © Copyright OHAUS Corporation

Click or scan for the web

