

Code: 22002073



TECHNICAL DATA SHEET

Performance

Heating power	600
Over Temperature Safety Device	420°C
Working temperature	5 ... 40
Heating Temperature Range	Ta ... 340
Working humidity	80 % RH
Stirring speed	50 ... 1500
Max stirring quantity	20 L
Display	LCD

Features

IP rating	42
Temperature indication	LCD display (°C resolution = 0,1°C)
Speed indication	LCD Display
Maximum length of magnetic anchor	80
Working plate - material	Ceramic-coated aluminum

Feeds

Power supply	220
Frequency	50 / 60

Dimensions and weight

Dimensions (WxDxH)	270x175x98
Weight [Kg]	2,2
Working plate - dimensions	135 diameter mm