

Code: 22003033



TECHNICAL DATA SHEET

Performance

Heating power	1150
Over Temperature Safety Device	580°C
Working temperature	5 ... 40
Heating Temperature Range	Ta ... 550 °C
Working humidity	80 % RH
Stirring speed	50 ... 1500
Max stirring quantity	20 L
Display	LCD

Features

IP rating	21
Maximum length of magnetic anchor	80
Working plate - material	Glass ceramic

Feeds

Power supply	100 - 240
Frequency	50 / 60

Dimensions and weight

Dimensions (WxDxH)	360 x 230 x 110
Weight [Kg]	4.3
Working plate - dimensions	184 x 184 mm
