

Code: 22003023



## TECHNICAL DATA SHEET

### Performance

Heating power	1150
Over Temperature Safety Device	580°C
Working temperature	5 ... 40
Heating Temperature Range	Ta ... 550 °C
Working humidity	80 % RH
Stirring speed	50 ... 1500
Max stirring quantity	20 L
Display	LCD

### Features

IP rating	21
Temperature indication	LCD display (°C resolution = 0,1°C)
Speed indication	LCD Display
Maximum length of magnetic anchor	80
Working plate - material	Glass ceramic

### Feeds

Power supply	220
Frequency	50 / 60

### Dimensions and weight

Dimensions (WxDxH)	360 x 230 x 110
Weight [Kg]	4,3
Working plate - dimensions	184 x 184 mm