



designed for scientists



R 1342 Propeller stirrer, 4-bladed

/// Data Sheet

Standard stirring element, 4-bladed. For drawing the material to be mixed from the bottom to the top. Local shearing forces. Generates axial flow in the vessel. Used at medium to high speeds.

Stirrer Ø: 50 mm

Shaft Ø: 8 mm

Shaft length: 350 mm



designed for scientists

Technical Data

Stirrer diameter [mm]	50
Shaft diameter [mm]	8
Shaft length [mm]	350
Speed max. [rpm]	2000
Material in contact with medium	AISI 316L
Weight [kg]	0.142