



designed for scientists



IKA Tip 10 µl box (960 pcs.)

/// Data Sheet

Universal fit pipette tips, 10 µl size, non-sterile.
Fits to all grey IKA pipettes.

Made of polypropylene and can be autoclaved.
Excellent volume accuracy due to precise tip orifice.
Minimal residual wetting due to controlled quality of raw materials used.
Guaranteed dust-free due to special ionisation treatment.



designed for scientists

96 pcs/box, 10 boxes included in scope of delivery.



designed for scientists