



designed for scientists



IKA Tip 200 µl box (960 pcs.)

/// Data Sheet

Universal fit pipette tips, 200 µl size, non-sterile.
Fits to all yellow IKA pipettes.

- Made of polypropylene and can be autoclaved.
- Excellent volume accuracy due to precise tip orifice.
- Minimal residual wetting due to controlled quality of raw materials used.
- Guaranteed dust-free due to special ionization treatment.



designed for scientists

** Recommended for use with IKA single channel pipettes **

96 pcs/box, 10 boxes included in scope of delivery.



designed for scientists