

M2-A Agitator ARGOLab magn.heating analogue

Code: 22002013



TECHNICAL DATA SHEET

Performance

Heating power	530
Over Temperature Safety Device	420
Working temperature	5 ... 40 °C
Heating Temperature Range	Ta ... 340
Working humidity	80 % RH
Stirring speed	0 ... 1500
Max stirring quantity	20 L

Features

IP rating	21
Temperature indication	Analogue
Speed indication	Analogue
Maximum length of magnetic anchor	80
Working plate - material	Stainless Steel

Feeds

Power supply	100 – 240
Frequency	50 / 60

Dimensions and weight

Dimensions (WxDxH)	280 x 160 x 100
Weight [Kg]	2.8
Working plate - dimensions	Ø135 mm