

Vacuum Pump VP-20

Code: 41203152



TECHNICAL DATA SHEET

Standards applied	
Regulations	2014/30/EU, 2014/35/EU, 2011/65/EU, 2006/42/EC, EN 60034-1:2010+AC:2010, EN IEC 61000-3-2:2019+A1:2021+A2:2024, EN 61000-3-3:2013+A2:2021+AC:2022, EN IEC 55014-1:2021 EN IEC 55014-2:2021
Warranty [years]	2 anni
Supplied as standard	
Supplied as standard	1 connecting pipe of 1 m
Performance	
Maximum Pump Flow Rate [l/min]	20 l/min
Ultimate vacuum	8 mbar
Speed (RPM)	1450 rpm
Wetted materia	PTFE, FFPM
Pump	2-stage pump
Dimensions and weight	
External dimension (WxDxH)	165 x 315 x 270 mm
Packaging dimension (WxDxH)	460 x 330 x 295 mm
Packaging weight (Kg)	12 kg