

Vacuum Controller VC-50

Code: 41203162



TECHNICAL DATA SHEET

Standards applied	
Warranty [years]	2 anni
Supplied as standard	
Supplied as standard	Kit of 2 connection pipes of 0,5 m each
Performance	
Control range	1 ... 1013 mbar
Display	Touchscreen
Control	Fixed value automatic control, preset program automatic control
Sensor	Ceramic
Dimensions and weight	
External dimension (WxDxH)	105 x 180 x 185 mm
Packaging dimension (WxDxH)	410 x 270 x 190 mm
Packaging weight (Kg)	3 kg