

C3-20 Chiller

Code: 41203112



TECHNICAL DATA SHEET

Standards applied

Regulations

2014/30/EU, 2014/35/EU, 2011/65/EU, 2006/42/EC, EN 61010-1:2010+A1:2019, EN 61010-2-010:2014, N 61010-2-011:2017, EN 61326-1:2013, EN 6100-3-2:2019, EN 6100-3-3:2013+A1:2019, UE regulation 517/2014 concerning F-GAS, Directive 2014/68 UE regarding PED

Warranty [years]

2 anni

Supplied as standard

Supplied as standard

Set of 2 x 1,5m pipes with insulating material

Performance

Temperature Range	-20 ... Ambient Temperature °C
Capacity	5 L
Temperature Resolution	0,1 °C
Volume	5 L
Display	LCD
Timer	99:59 or continuous mode hh:mm
Security Class (DIN 12880)	2.0
Circulation Type	Forced circulation
Pump	Yes
Maximum Pump Flow Rate [l/min]	6 l/min
Max head	9 meters
Power	550 W
Overtemperature Safety System	Electronic
Refrigerant	R 290
External Circulation	Yes
Refrigeration Capacity at 10°C	530 W

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

Power Supply	230 V
--------------	-------

Dimensions and weight

Weight	32 kg
External dimension (WxDxH)	275 x 520 x 610 mm
Packaging dimension (WxDxH)	675 x 430 x 745 mm
Packaging weight (Kg)	50 kg