

WITOPET professional

The professional way for your pipetting applications

WITOPET professional single channel is the classical microliter pipette, made for the daily use in the laboratory: accurate, robust and handy. The seven models available cover a volume range from 0,2 μ l to 5 ml.

Push-button

The push button is color-coded to avoid confusion and softly gliding for comfortable pipetting.

Adjustment screw

The color coded adjustment screw ensures easy and accurate volume setting.



Color-Code

The push button, the adjustment thread and the tip ejector button are color coded and help to distinguish all pipettes at a single glance.



Piston

The optimized piston allows fine adjustment of the pipetting forces, thus making the WITOPET professional perfect for precise laboratory applications like e.g. gel loading.



Tip ejector

The tip ejector fits well into the body of the pipette to avoid slipping over time. It is made of stainless steel for a long lifetime.



Tip cone

The tip cone is suitable for all established brands of pipette tips.



Single-Channel Microliter Pipette

Your advantages at a glance:

- high quality pipette, made in Germany
- economical priced
- perfect for precise laboratory applications
- 7 color-coded single channel models
- high accuracy and precision according to DIN EN ISO 8655
- ergonomic designed and low weighted
- optimized spring force for comfortable pipetting
- robust and stable construction
- tip ejector and piston made of stainless steel
- suitable for all established brands of pipette tips
- completely autoclavable at 121°C
- individual certificate and serial number
- DE-M marked according to the German calibration law, valid from January 2015

WITOPET professional

Volumen µl	Color	Division µl	Volume tested µl	Accuracy ± µl	Coefficient of variation ≤ µl	Order number
0.2 - 2.0	red	0.1	0.2	0.08	0.04	5 402 002
			0.5	0.08	0.04	
			1	0.08	0.04	
			2	0.08	0.04	
0.5 - 10	red	0.1	1	0.12	0.08	5 402 010
			5	0.12	0.08	
			10	0.12	0.08	
2 - 20	yellow	0.1	2	0.20	0.10	5 402 020
			10	0.20	0.10	
			20	0.20	0.10	
10 - 100	yellow	1.0	10	0.80	0.30	5 402 100
			50	0.80	0.30	
			100	0.80	0.30	
20 - 200	yellow	1.0	20	1.60	0.60	5 402 200
			100	1.60	0.60	
			200	1.60	0.60	
100 - 1000	blue	10	100	8	3	5 402 901
			500	8	3	
			1000	8	3	
1000 - 5000	black	10	1000	40	15	5 402 950
			2500	40	15	
			5000	40	15	
Pipette support for Witopet professional pipettes						5 406 107

