

# aeroject® ultra Pipet Tips

free of DNA, DNase,  
pyrogens and ATP








## aeroject® ultra tips – Pipet Tips with Aerosol Barrier

CE IVD

Pipet tips with special filters, to effectively prevent the pipettor shaft from contamination with aerosols and liquids to be sucked into the pipettor.

- sterilized
- free of RNase-, DNase-, DNA-, pyrogens and ATP
- free of human DNA and PCR inhibitors
- pipet tips and aerosol filters made of ultra-hydrophobic HD-PE
- minimized retention due to special surface treatment
- perfect and tight fit on all major pipettor brands
- aeroject® ultra 17 70 010 and 17 71 000xl with extended length for operation with narrow vessels and reduced risk of cross-contamination
- packed in plastic racks with hinged lids

Typ	Volume Range $\mu$ l	Length mm	Packaging	Quantity per Pack	Order No.
Pipet Tip aeroject® ultra tips xl, for all commonly used pipets, ratiopetta®	0,1-10	45,20	Racks 10 x 96	960	17 70 010
					
Pipet Tip aeroject® ultra tips, for all commonly used pipets, ratiopetta®	0,1-20	50,80	Racks 10 x 96	960	17 70 020
					
Pipet Tip aeroject® ultra tips, for all commonly used pipets, ratiopetta®	1-100	50,80	Racks 10 x 96	960	17 70 100
					
Pipet Tip aeroject® ultra tips, for all commonly used pipets, ratiopetta®	1-200	50,80	Racks 10 x 96	960	17 70 200
					
Pipet Tip aeroject® ultra tips, for all commonly used pipets, ratiopetta®	100-1.000	70,10	Racks 10 x 60	600	17 71 000
					
Pipet Tip aeroject® ultra tips, for all commonly used pipets, ratiopetta®	100-1.000	86,00	Racks 10 x 96	960	17 71 000xl
